



For immediate release:

February 11, 2015

Sean Gilmore

Distek, Inc.

121 North Center Drive

North Brunswick, NJ 08902

T: 732 422 7585

E: info@distekinc.com

Distek, Inc. Releases Mini HeightChek™ Tool for Small Volume Dissolution

North Brunswick, NJ – Distek, Inc., an industry leader in pharmaceutical laboratory testing instruments, accessories and validation services, releases the Mini HeightChek Tool for Small Volume Dissolution.

The Mini HeightChek Tool is a fixed 10mm gauge that accurately sets the paddle or basket heights prior to using the Distek small volume dissolution accessory.

“As we’ve see an increase in utilization of small volume (<200 ml) dissolution testing in the marketplace we are pleased to be able to provide our customers with a quick and simple tool for setting the operating heights,” says Jeff Seely, Vice President of Sales & Business Development, Distek, Inc.

For additional information contact Distek’s Corporate Office at +1 732 422 7585, by email at info@distekinc.com or visit www.distekinc.com.

ABOUT DISTEK, INC.

Distek, Inc., headquartered in North Brunswick, NJ, is a leading manufacturer of pharmaceutical laboratory test equipment, specializing in dissolution products and services. Distek also provides solutions for UV fiber optics, media preparation, physical testing, disintegration, and validation services.



**Now Available
Mini HeightChek Tool
for Small Volume
Dissolution**